

TEKIA

Injector Recommendations

Validation #	Injector	Cartridge	Tekia Lens Models
2009-001	Medicel VISCOJET	Medicel Microglide 2.8 mm	500
2007-004	Medicel Viscojet	Medicel Microglide 2.8 mm	811, 820, 872 MeyeLens 260,460 Ophthalmo Pro AC 3000 HD
2016-007	Medicel Accuject Pro screw 2.2S-BL	Medicel Microglide 2.2 mm	900
2006-007	Medicel NAVIJECT	Medicel NAVIGLIDE 2.9 - 3 P	617 618 MeyeLens 360