

Indications: Posterior chamber/cataract surgery


### Additional Specifications

Not available in the United States

Optical hydrophobic acrylic with UV blocker  
 Diopter range: 10.0 to 30.0 @ 0.5 increments  
 Refractive Index: 1.555 @ 35°C

Light transmission (380 – 800 nm): > 96%  
 UV cut-off < 10%: 379 nm  
 UV transmission (300 – 380 nm): 0.5%

Biconvex aspheric design  
 360° optical square edge  
 Incision size: 2.2 mm



A-Constant	Ultrasound A-Scan	Optical Biometry
SRK/T Formula	117.8	118.8
SRK II Formula	118.1	119.1

# TEKIA

